

Westlab's Dissection Range

Codes 663-925 to 663-938



Westlab's Dissection Range

Westlab's Dissection Products are manufactured from a Surgical Grade Stainless Steel Material to avoid corrosion or rust. All the conventional sterilisation methods will not cause damage to these instruments.

The range is suitable for use in the classroom by teachers and students alike.



Product Features

- Material: Surgical Grade Stainless Steel
- Non-Sterile, Reusable
- Made In Pakistan

PRODUCT CODES

663-925 to
663-938

Australia

1800 358 101
sales.au@westlab.com
www.westlab.com.au

Canada

1 877 822 1455
sales.ca@westlab.com
www.westlab.com

New Zealand

+64 9 553 3677
sales.nz@westlab.com
www.westlab.co.nz

United States

646 568 5391
sales.us@westlab.com
www.westlab.com

Dissecting Kits

Product Specifications

663-925 Dissecting Kit, 7 Piece

Kit Items	Qty
Dressing Thumb Forceps, Straight, 130mm	1
Dressing Thumb Forceps, Straight, 115mm	1
Dressing Scissors Sharp/Blunt, Straight, 130mm	1
Iris Scissors, Straight, 115mm	1
Scalpel Handle #3	1
Scalpel Blade #10	1
Wood Needle	1
Probe	1
Black Pouch for Dissection Set	1

663-926 Dissecting Kit, 14 Piece

Kit Items	Qty
Dressing Thumb Forceps, Straight, 130mm	1
Dressing Thumb Forceps, Straight, 115mm	1
Dressing Scissors Sharp/Blunt, Straight, 130mm	1
Iris Scissors, Straight, 115mm	1
Dissecting Needle	1
Scalpel Handle #3	1
Scalpel Blade #10	1
Arrow Point Needle	1
Magnifying Glass Mirror	1
Dissecting Hook	1
Blow Pipe	1
Spatula	1
Razor	1
Black Pouch for Dissection Set	1



PRODUCT CODES


663-925
663-926



Dissecting Kit, 7 Piece



Dissecting Kit, 14 Piece



Dissecting Forceps

Product Specifications

Codes	Description	Length	Pack Size
663-931	Forceps, Artery Straight, Stainless Steel	140mm	10
663-932	Forceps, Artery Straight, Stainless Steel	200mm	10
663-935	Forceps, Thumb Straight, Blunt Ends, Stainless Steel	130mm	10
663-936	Forceps, Thumb Straight, Blunt Ends, Stainless Steel	200mm	10
663-937	Forceps, Microscopic, Pointed Ends, Stainless Steel	105mm	10
663-938	Forceps, Thumb Curved, Pointed Ends, Stainless Steel	130mm	10

PRODUCT CODES

663-931
 663-932
 663-935
 663-936
 663-937
 663-938

Australia

1800 358 101
sales.au@westlab.com
www.westlab.com.au

Canada

1 877 822 1455
sales.ca@westlab.com
www.westlab.com

New Zealand

+64 9 553 3677
sales.nz@westlab.com
www.westlab.co.nz

United States

646 568 5391
sales.us@westlab.com
www.westlab.com




Dissecting Scissors

Product Specifications

Codes	Description	Length	Pack Size
663-933	Scissors, Sharp/Blunt, Stainless Steel	130mm	10
663-934	Scissors, Sharp/Sharp, Stainless Steel	130mm	10

PRODUCT CODES

663-933
663-934



Miscellaneous Dissecting Accessories

Product Specifications

Codes	Description	Length	Pack Size
663-927	Blow Pipe, Stainless Steel	200mm	10
663-928	Dissecting Needles, Stainless Steel, Pair of Sharp (1) & Blunt (1)	127mm	10
663-929	Blades Scalpel, Handle No. 3	126mm	10
663-930	Blades Scalpel, Handle No. 4	137mm	10

PRODUCT CODES

663-927
 663-928
 663-929
 663-930

Australia

1800 358 101
sales.au@westlab.com
www.westlab.com.au

Canada

1 877 822 1455
sales.ca@westlab.com
www.westlab.com

New Zealand

+64 9 553 3677
sales.nz@westlab.com
www.westlab.co.nz

United States

646 568 5391
sales.us@westlab.com
www.westlab.com