




Westlab Precision Balances

Codes 665-080, 665-081, 665-082, 665-083, 665-084



Westlab Precision Balances

Designed with high precision load cells and solid structures, Westlab's Precision Balances guarantee reliability and stability. The stainless-steel pans are corrosion resistant and easy to clean.

The integrated rechargeable battery allows the user to move from lab-to-lab or use the unit in an area with no access to a power source. In a teaching environment, science educators can charge the balances before a lesson to avoid the power cords being tangled. Additionally, the low or charged battery indicator displayed on the screen ensures that sample weighing is not interrupted.

With a backlit LCD screen, the balances feature a simple four key menu for easy navigation.

Product Features

- LCD Display, Backlit
- High Precision Load Cell
- Integrated Rechargeable Battery
- Selectable Measuring Units
- Overload Display Indicator
- Piece Counting Function
- Height Adjustable Feet
- Intelligent Calibration
- Full Capacity Taring

PRODUCT CODES

665-080
665-081
665-082
665-083
665-084

Australia

1800 358 101
sales.au@westlab.com
www.westlab.com.au

Canada


1 877 822 1455
sales.ca@westlab.com
www.westlab.com

New Zealand

+64 9 553 3677
sales.nz@westlab.com
www.westlab.co.nz

United States

646 568 5391
sales.us@westlab.com
www.westlab.com



Precision Balance 300g x 0.001g

Product Specifications

Display	LCD, backlit
Load Rating	300g
Accuracy	0.001g
Tare Range	0 – 300g
Repeatability	1mg
Linearity	± 3mg
Overload Protection	Yes
Feet	Height adjustable
Power Supply	220 – 240V or rechargeable battery
Pan Diameter	110mm
Unit Dimensions	280mm L x 180mm W x 180mm H

PRODUCT CODES

665-080

Australia

1800 358 101
sales.au@westlab.com
www.westlab.com.au

Canada

1 877 822 1455
sales.ca@westlab.com
www.westlab.com

New Zealand

+64 9 553 3677
sales.nz@westlab.com
www.westlab.co.nz

United States

646 568 5391
sales.us@westlab.com
www.westlab.com



Precision Balance 600g x 0.1g

Product Specifications

Display	LCD, backlit
Load Rating	600g
Accuracy	0.1g
Tare Range	0 – 600g
Repeatability	0.1g
Linearity	± 0.2g
Overload Protection	Yes
Feet	Height adjustable
Power Supply	220 – 240V or rechargeable battery
Pan Dimensions	143mm x 192mm
Unit Dimensions	220mm L x 220mm W x 70mm H

PRODUCT CODES

665-081

Australia

1800 358 101
sales.au@westlab.com
www.westlab.com.au

Canada

1 877 822 1455
sales.ca@westlab.com
www.westlab.com

New Zealand

+64 9 553 3677
sales.nz@westlab.com
www.westlab.co.nz

United States

646 568 5391
sales.us@westlab.com
www.westlab.com



Precision Balance 3000g x 0.01g

Product Specifications

Display	LCD, backlit
Load Rating	3000g
Accuracy	0.01g
Tare Range	0 - 3000g
Repeatability	0.1g
Linearity	± 0.03g
Overload Protection	Yes
Feet	Height adjustable
Power Supply	220 – 240V or rechargeable battery
Pan Dimensions	200mm x 140mm
Unit Dimensions	220mm L x 220mm W x 70mm H

PRODUCT CODES

665-082

Australia

1800 358 101
sales.au@westlab.com
www.westlab.com.au

Canada

1 877 822 1455
sales.ca@westlab.com
www.westlab.com

New Zealand

+64 9 553 3677
sales.nz@westlab.com
www.westlab.co.nz

United States

646 568 5391
sales.us@westlab.com
www.westlab.com



Precision Balance 6000g x 0.1g

Product Specifications

Display	LCD, backlit
Load Rating	6000g
Accuracy	0.1g
Tare Range	0 – 6000g
Repeatability	0.1g
Linearity	± 0.3g
Overload Protection	Yes
Feet	Height adjustable
Power Supply	220 – 240V or rechargeable battery
Pan Dimensions	143mm x 192mm
Unit Dimensions	220mm L x 220mm W x 70mm H

PRODUCT CODES

665-083

Australia

1800 358 101
sales.au@westlab.com
www.westlab.com.au

Canada

1 877 822 1455
sales.ca@westlab.com
www.westlab.com

New Zealand

+64 9 553 3677
sales.nz@westlab.com
www.westlab.co.nz

United States

646 568 5391
sales.us@westlab.com
www.westlab.com



Precision Balance 10kg x 0.1g

Product Specifications

Display	LCD, backlit
Load Rating	10,000g
Accuracy	0.1g
Tare Range	0 – 10,000g
Repeatability	± 0.3g
Linearity	± 0.4g
Overload Protection	Yes
Feet	Height adjustable
Power Supply	220 – 240V or rechargeable battery
Pan Dimensions	170mm x 170mm
Unit Dimensions	280mm L x 180mm W x 80mm H

PRODUCT CODES

665-084

Australia

1800 358 101
sales.au@westlab.com
www.westlab.com.au

Canada

1 877 822 1455
sales.ca@westlab.com
www.westlab.com

New Zealand

+64 9 553 3677
sales.nz@westlab.com
www.westlab.co.nz

United States

646 568 5391
sales.us@westlab.com
www.westlab.com