

# Statement of Compliance

## Westlab Safety Cabinet Compliance

Products: Indoor safety cabinets of various classes.

### Statement

To whom it may concern, Westlab manufactures in Australia a range of safety cabinets for storage of various classes of dangerous goods to exceed the requirements of the applicable Australian Standard for each class to be stored.

All safety cabinets manufactured comply in functionality, construction, performance and appearance to the following Australian Standards:

- Class 3 Flammable Liquids Storage: AS1940-2004
- Class 5.1 Oxidising Agents Storage: AS4326-2008
- Class 5.2 Organic Peroxides Storage: AS2714-2008
- Class 6 Toxic Substances Storage: AS/NZS4452:1997
- Class 8 Corrosive Substances Storage: AS3780-2008

All double door cabinets (steel construction) feature the exclusive SAFE-T-CLOSE sequential closing system which guarantees effective, simple and trouble free automatic door closing.

### Westlab codes covered by these standards

663-100	30L	4-in-1-Safety Cabinet
663-101	60L	4-in-1-Safety Cabinet
663-102	100L	4-in-1-Safety Cabinet
663-103	160L	4-in-1-Safety Cabinet
663-104	250L	4-in-1-Safety Cabinet

### Authorised by

**Ted Fowler**  
Technical Director

**Date:** 07/28/2020