



Westlab Safety Glasses


Codes 663-820, 663-821, 663-822



Safety Glasses, Discovery

Westlab's Discovery Safety Glasses are the most common type used in Australian laboratories. With non-adjustable arms and good eye coverage, these are an affordable quality pair of goggles.

These glasses are ideal for humid conditions due to their anti-fog coated lens. Additionally, the hard coated lenses offer anti-scratch protection. The lens also filters out 99.9% of UV radiation.



Product Features

- Sleek & Ultra Lightweight Wraparound Lens
- Anti-Fog Coated Lens, Ideal for Humid Conditions
- AS/NZS 1337.1:2010 Certified, BMP 557823
- Hard Coated Lens for Anti-Scratch Protection
- 99.9% UV Protection

Product Specifications

Colour	Clear
Purchase Quantity	Each
Box Quantity	12
Carton Quantity	300
UV Protection	99.9%
Certification	AS/NZS 1337.1:2010 BMP 557823

PRODUCT CODE

663-820

Australia

1800 358 101
sales.au@westlab.com
www.westlab.com.au

Canada

1 877 822 1455
sales.ca@westlab.com
www.westlab.com

New Zealand

+64 9 553 3677
sales.nz@westlab.com
www.westlab.co.nz

United States

646 568 5391
sales.us@westlab.com
www.westlab.com




Safety Glasses, Research

Westlab's Research Safety Glasses are designed to be worn over prescription glasses. These glasses are ideal for humid conditions due to their anti-fog coated lens.

Additionally, the hard coated lenses offer anti-scratch protection. The lens also filters out 99.9% of UV radiation.

Suitable for site visitors, school children, laboratories, and prescription glasses wearers.



Product Features

- Suitable for Prescription Glasses
- Integrated Side Shields, Vented
- Anti-Fog Coated Lens, Ideal for Humid Conditions
- AS/NZS 1337.1:2010 Certified, BMP 557823
- 99.9% UV Protection

Product Specifications

Colour	Clear
Purchase Quantity	Each
Box Quantity	12
Carton Quantity	300
UV Protection	99.9%
Certification	AS/NZS 1337.1:2010 BMP 557823

PRODUCT CODE

663-821

Australia

1800 358 101
sales.au@westlab.com
www.westlab.com.au

Canada

1 877 822 1455
sales.ca@westlab.com
www.westlab.com

New Zealand

+64 9 553 3677
sales.nz@westlab.com
www.westlab.co.nz

United States


646 568 5391
sales.us@westlab.com
www.westlab.com



Safety Glasses, Lab-Protect

Westlab's Lab-Protect Safety Glasses have a classic lightweight design and length adjustable side arms. These glasses are ideal for humid conditions due to their anti-fog coated lens.

Additionally, the hard coated lenses offer anti-scratch protection. The lens also filters out 99.9% of UV radiation for outdoor use.



Product Features

- Classic & Popular Lightweight Design
- Anti-Fog Coated Lens, Ideal for Humid Conditions
- Hard Coated Lens for Anti-Scratch Protection
- AS/NZS 1337.1:1992 Certified
- Length Adjustable Side Arms
- 99.9% UV Protection

Product Specifications

Colour	Clear
Purchase Quantity	Each
Box Quantity	12
Carton Quantity	300
UV Protection	99.9%
Certification	AS/NZS 1337.1:1992

PRODUCT CODE

663-822

Australia

1800 358 101
sales.au@westlab.com
www.westlab.com.au

Canada

1 877 822 1455
sales.ca@westlab.com
www.westlab.com

New Zealand

+64 9 553 3677
sales.nz@westlab.com
www.westlab.co.nz

United States

646 568 5391
sales.us@westlab.com
www.westlab.com