


Dip & Read Tests

For the Food Industry



New: Precise readings with QUANTOFIX® Relax



Rapid Tests for the Food Industry

- Fast & easy
- Reliable results
- Efficient quality control

MACHERY-NAGEL

www.mn-net.com









MN Dip & Read Tests for the Food Industry

High quality and consumer safety is a fundamental aspect in food production. MN offers a wide range of test strips that can be used to verify constant high quality of products and safety of equipment and production facilities. MN test strip solutions are fast, easy to use and reliable, making products and operations safe and efficient at the same time.



Product information

Test	Product	Application	Gradation	REF
	QUANTOFIX® Peroxide 25	<ul style="list-style-type: none"> Reliable detection of peroxide residues down to 0.5 mg/L Can be stored at room temperature Accurate results in 15 seconds 	0 · 0.5 · 2 · 5 · 10 · 25 mg/L H ₂ O ₂	913 19
	QUANTOFIX® Sulfite	<ul style="list-style-type: none"> Sulphite in food can cause severe allergic reactions Tests for sulphite in food treated with sulphur based preservatives Reliable results in 20 seconds 	0 · 10 · 25 · 50 · 100 · 250 · 500 · 1000 mg/L SO ₃ ²⁻	913 06
	QUANTOFIX® Peracetic Acid 2000	<ul style="list-style-type: none"> Specific to peracetic acid No peroxide interference during reaction time Detects disinfectant potency levels up to 2000mg/L 	0 · 500 · 1000 · 1500 · 2000 mg/L peracetic acid	913 42
	pH-Fix 2.0-9.0	<ul style="list-style-type: none"> Always ready-to-use Brilliant color chart for precise results Reliable results in seconds 	pH 2.0 · 2.5 · 3.0 · 3.5 · 4.0 · 4.5 · 5.0 · 5.5 · 6.0 · 6.5 · 7.0 · 7.5 · 8.0 · 8.5 · 9.0	921 18
	QUANTOFIX® Ascorbic Acid (Vitamin C)	<ul style="list-style-type: none"> Vitamin C is important and desirable in food Also used as stabilizer and reducing agent Reliable results in 45 seconds 	0 · 50 · 100 · 200 · 300 · 500 · 1000 · 2000 mg/L vitamin C	913 14
	QUANTOFIX® Nitrate/Nitrite	<ul style="list-style-type: none"> Prevalent in many foods High doses have adverse health effects Reliable results in 60 seconds 	0 · 10 · 25 · 50 · 100 · 250 · 500 mg/L NO ₃ ⁻ 0 · 1 · 5 · 10 · 20 · 40 · 80 mg/L NO ₂ ⁻	913 13
Others				
	Peroxtesmo KO	Quick and easy detection of peroxidase in food		906 06
	PEHANON	Easy pH determination in colored liquids		904 ...
	QUANTOFIX® Cyanide	Reliable cyanide detection, especially in fruit brandy production		913 18
	INDIPRO	Detection of protein residue on equipment and working surfaces		907 65

QUANTOFIX® Relax

REF 913 46

The instrument is for the automatic reading of QUANTOFIX® test strips. Just place the strip on the strip adaptor and the Relax will give a precise and reliable result within seconds. The results are printed, saved internally and can be transmitted to a PC. It is ideal for QC managers and decision makers in production.

QUANTOFIX® Relax – for easy, reliable QC procedures

- highest precision
- objective results
- easy result management



These tests are suitable for detection with QUANTOFIX® Relax.



www.mn-net.com

MACHEREY-NAGEL



MACHEREY-NAGEL GmbH & Co. KG · Neumann-Neander-Str. 6–8 · 52355 Düren · Germany
Germany and international:
 Tel.: +49 24 21 969-0
 Fax: +49 24 21 969-199
 E-mail: info@mn-net.com

Switzerland:
MACHEREY-NAGEL AG
 Tel.: +41 62 388 55 00
 Fax: +41 62 388 55 05
 E-mail: sales-ch@mn-net.com

France:
MACHEREY-NAGEL EURL
 Tel.: +33 388 68 22 68
 Fax: +33 388 51 76 88
 E-mail: sales-fr@mn-net.com

USA:
MACHEREY-NAGEL Inc.
 Tel.: +1 484 821 0984
 Fax: +1 484 821 1272
 E-mail: sales-us@mn-net.com

